

TEXAS A&M UNIVERSITY

PHYSICS LECTURE DEMONSTRATIONS

Colleagues,

The Physics Demonstration Program has been reorganized and updated for implementation in the new George P. and Cynthia Woods Mitchell Institute for Fundamental Physics and Astronomy and the George P. Mitchell '40 Physics Building. [Click here to view Optics Demonstrations.](#)

The labeling system of the demonstrations has also been changed. The new labeling will incorporate the Demonstration Classification System (DCS) suggested by PIRA, the Physics Instructional Resource Association. The goal of the PIRA Demonstration Classification Scheme is to create a logically organized and universally inclusive taxonomy giving a unique number to every lecture demonstration.

External website references have been added to the demonstration notes when possible. These references have information (historical and practical) to supplement the notes listed in the demo page.

If you are interested in using one of these demonstrations in your classroom, you may send me an email with your class information and the demonstration number and name. Please allow at least two class days notice for the demonstration requests.

Santos Ramirez

Coordinator

-- Physics Lecture Demonstration Program

-- 200 Level Physics Teaching Laboratory Program

Department of Physics and Astronomy

Texas A&M University

College Station, TX 77843-4242

U.S.A.

Phone: (979) 458-7918

Email: Ramirez@physics.tamu.edu

Optics

Geometric Optics

6A20.35 MIRAGE (REFLECTION FROM CURVED SURFACES)

6A44.40 LIGHT THROUGH A THICK FIBER OPTIC

6A60.10 WHITE BOARD OPTICS (CONCAVE AND CONVEX LENSES)

6A60.30 IMAGE FORMATION WITH LENSES

Interference

6D10.10 WHITE LIGHT AND CORNELL SLIDES

6D20.15 DIFFRACTION GRATING / LASER

6D20.20 DIFFRACTION GRATING / SPECTRUM OF DISCHARGE TUBES

Color

6F30.10 DISPERSION OF LIGHT BY A PRISM


Polarization

6H10.10 POLAROIDS ON AN OVERHEAD

6H10.20 MICROWAVE POLARIZATION

6H50.00 DEPOLARIZATION BY SCATTERING USING WAXPAPER

6A20.35 MIRAGE (REFLECTION FROM CURVED SURFACES)




Apparatus: Mirage mirrors and a small pig

Notes: Inside the mirage are two concave mirrors opposite each other. The image (a small pig) is placed inside. This soft figure is used to prevent scratching the mirror. The pig appears to be floating on top of the mirage fixture. A web cam could be used to make this demo more visible.

References: <http://spie.org/etop/2007/etop07methodsV.pdf>

[Top of the Document](#)

6A44.40 LIGHT THROUGH A THICK FIBER OPTIC




Apparatus: Clear curved lucite rod on stand
Laser pointer

Notes: The total internal refraction of the lucite rod may be shown.

DO NOT SHINE THE LASER LIGHT TOWARD THE CLASS.

[Top of the Document](#)

6A60.10 WHITE BOARD OPTICS (CONCAVE AND CONVEX LENSES)



Apparatus: Ray tracing apparatus (light and adaptor)
Lenses and prisms mounted on apparatus

Notes: Attach the lens or prisms on the metal screen (with the mounted magnets). The light source may be adjusted for a single light ray or multiple light rays.

[Top of the Document](#)

6A60.30 IMAGE FORMATION WITH LENSES



Apparatus: Optical bench
incandescent lamp
converging lens
diverging lens
frosted plate

Notes: Plug the lamp into the wall outlet and adjust the distance between the lens and the lamp to form an image on the frosted glass plate. The frosted plate should be toward the students. Is the image upside down? Focal length and magnification of a lens may be shown.

[Top of the Document](#)

6D10.10 WHITE LIGHT AND CORNELL SLIDES



Apparatus: Light with a vertical filament
Box of Cornell slides
Cornell slide sheet (that shows the different patterns on the slide)

Notes: The students can look through the different slits in the Cornell slide to see diffraction and interference patterns. (the slide has single slits, double slits and 4 different diffraction gratings)

DO NOT SHINE THE LASER LIGHT TOWARD THE CLASS.

[Top of the Document](#)

6D20.15 DIFFRACTION GRATING / LASER



Apparatus: Cornell slide on stand
Laser on lab jack
Cornell slide sheet (that shows the different patterns on the slide)

Notes: Shine the laser at the diffraction grating on the Cornell slide. These are located on the middle row of the Cornell slide. Raise or lower the lab jack to select the other 4 different gratings on the slide.

References: http://www.chemie.de/lexikon/e/Helium-neon_laser/

DO NOT SHINE THE LASER LIGHT TOWARD THE CLASS.

[Top of the Document](#)

6D20.20 DIFFRACTION GRATING / SPECTRUM OF DISCHARGE TUBES



Apparatus: 2 gas discharge tubes in metal boxes
White light with vertical filament
Box of diffraction gratings

Notes: Stand the discharge tube boxes and white light on the desk top or demo table.
Students may use the diffraction gratings to observe the spectrum of the different gases and of the white light.

[Top of the Document](#)

6F30.10 DISPERSION OF LIGHT BY A PRISM



Apparatus: Prism on a mount
Overhead projector with two card board sheets

Notes: Adjust the card board sheets on the overhead projector for a slit of about 2 mm. Shine the light on the prism and adjust the prism such that a spectrum of colors is displayed on the wall or screen. The height of the prism must be in line with the head of the projector.

[Top of the Document](#)

6H10.10 POLARIZING SHEETS ON AN OVERHEAD PROJECTOR



Apparatus: Polaroid sheets (2)
Plastic protractor
Transparency sheet with “magic” tape
Overhead projector

Notes: Place the polaroid sheets on the overhead projector and rotate one of the sheets. Note the polarization effect. Put the transparency sheet (with the magic tape) between the polaroid sheets. Rotate one of the polaroid sheets (to show how the light through the magic tape changes). Put the protractor between the polaroid sheets and slightly flex the protractor. This will show the stresses in the plastic (but be careful not to break the plastic).

[Top of the Document](#)

6H10.20 MICROWAVE POLARIZATION




Apparatus: Microwave transmitter and receiver assembly
Metal grid plate
Aluminum plate
Wood-fiber plate

Notes: On the transmitter unit, set the transmitter klystron voltage control to mid-range and the oscillator switch to internal position. On the receiver unit, adjust the receiver gain control to mid-range. The horns of the two units should be aimed at each other and about 12 inches apart. Turn on the two units and adjust the speaker volume control to a suitable level. Place the different plates in front of the horns to show polarization and reflection.

[Top of the Document](#)

6H50.00 DEPOLARIZATION BY SCATTERING USING WAXPAPER



Apparatus: Polaroid sheets (2)
Sheet of waxpaper
Overhead projector

Notes: Place the polaroid sheets on the overhead projector and place the sheet of waxpaper between them. Rotate the polaroid sheets and show how the waxpaper dispersion will not allow the light to get dark.

[Top of the Document](#)